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

TaqMan

Polymerization


Step 1: A reporter (R) and a quencher (Q) are attached to the 5' and 3' ends of a TaqMan probe.

Strand Displacement


Step 1 continued: when both dyes are attached to the probe, reporter dye emission is quenched.

Cleavage

Step 2: During each extension cycle, the AmpliTaq Gold[®] DNA polymerase cleaves the reporter dye from the probe.

Polymerization Completed

Step 3: After being separated from the quencher, the reporter dye emits its characteristic fluorescence.


SYBR Green


Step 1: Reaction setup
The SYBR[®] Green I dye fluoresces when bound to double-stranded DNA.

Step 2: Denaturation
When the DNA is denatured, the SYBR[®] Green I dye is released and the fluorescence is drastically reduced.

Step 3: Polymerization
During extension, primers anneal and PCR product is generated.

Step 4: Polymerization completed
SYBR[®] Green I dye binds to the double-stranded product, resulting in a net increase in fluorescence detected by the instrument.

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
If the efficiencies of the two systems are...	Then...
< 0.1	<ul style="list-style-type: none"> You can use the $\Delta\Delta C_T$ calculation for the relative quantification of target without running standard curves on the same plate.
> 0.1	<ul style="list-style-type: none"> Perform the validation experiment over a larger dynamic range (5 to 6 logs) or design and synthesize new primers to improve efficiency. Perform the quantification using the standard curve method.

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
Figure 2: Inhibition Plot Sample #1 demonstrating no inhibition


Figure 3: Inhibition Plot Sample #2 demonstrating presence of inhibitor

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RefGenes is a new type of tool that allows the identification of the genes having the highest stability of expression across a chosen set of conditions.

RefGenes selects candidate genes from the entire genome. The primary objective is to provide scientists performing qRT-PCR experiments an objective choice of reference genes specific for their experimental context.

The screenshot shows the GENEVESTIGATOR RefGenes tool interface. The interface is divided into several sections:

- SINGLE EXPERIMENT ANALYSIS**: Includes 'SINGLE EXPERIMENT TOOLS' with 'Samples' and 'Differential Expression'.
- COMPENDIUM-WIDE ANALYSIS**: Includes 'CONDITION SEARCH TOOLS' with 'Anatomy', 'Cell Lines', 'Cancers', 'Perturbations', and 'Development'.
- GENE SEARCH TOOLS**: Includes 'Anatomy', 'Cell Lines', 'Cancers', 'Perturbations', 'Development', 'RefGenes' (highlighted with a red circle), and 'Ortholog Search'.
- SIMILARITY SEARCH TOOLS**: Includes five tools for similarity analysis.

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

Times and Temperatures (Two-step RT-PCR)				
1) RT Step	HOLD	HOLD	* For reference only. RT is complete at this point.	
	10 min @ 25 °C	120 min @ 37 °C		
2) PCR Step	Initial Steps		PCR (Each of 40 cycles)	
	AmpErase® UNG Activation	AmpliTaq Gold® DNA Polymerase Activation	Melt	Anneal/Extend
	HOLD	HOLD	CYCLE	
	2 min @ 50 °C	10 min @ 95 °C	15 sec @ 95 °C	1 min @ 60 °C

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Primers	Selection Guidelines
Oligo d(T) ₁₆	<ul style="list-style-type: none"> • Use to reverse transcribe only eukaryotic mRNAs and retroviruses with poly-A tails • Avoid long mRNA transcripts or amplicons greater than 2 kilobases upstream
Random hexamers	<ul style="list-style-type: none"> • Try first for use with long reverse transcripts or reverse transcripts containing hairpin loops • Use to transcribe all RNA (rRNA, mRNA, and tRNA)
Sequence-specific reverse primers	<ul style="list-style-type: none"> • Use to reverse transcribe RNA-containing complementary sequences only

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Primer and probe design guidelines for quantitative assays

Probe Guidelines	Primer Guidelines
Select the probe first and design the primers as close as possible to the probe without overlapping it (amplicons of 50 to 150 basepairs are strongly recommended).	
Keep the G/C content in the 30 to 80% range.	
Avoid runs of an identical nucleotide. This is especially true for guanine, where runs of four or more Gs should be avoided.	
When using Primer Express software, the T _m should be 68 to 70 °C.	When using Primer Express software, the T _m should be 58 to 60 °C.
No G on the 5' end.	The five nucleotides at the 3' end should have no more than two G and/or C bases.
Select the strand that gives the probe more C than G bases.	
Make TaqMan MGB probes as short as possible, without being shorter than 13 nucleotides.	

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Recommended primer concentrations for DNA and cDNA quantification assays

Chemistry	Concentrations (nM)	
	Forward Primer	Reverse Primer
TaqMan probe	900	900
SYBR Green I dye	50	50

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Plot group key:

A: Combinations that contain at least 300 nM of forward and reverse primer

B: Combinations that contain at least 150 nM of forward and reverse primer

C: Combinations that contain at least 50 nM of forward and reverse primer

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Cebadores de la RT:

Primers	Selection Guidelines
Oligo d(T) ₁₆	<ul style="list-style-type: none">• Use to reverse transcribe only eukaryotic mRNAs and retroviruses with poly-A tails• Avoid long mRNA transcripts or amplicons greater than 2 kilobases upstream
Random hexamers	<ul style="list-style-type: none">• Try first for use with long reverse transcripts or reverse transcripts containing hairpin loops• Use to transcribe all RNA (rRNA, mRNA, and tRNA)
Sequence-specific reverse primers	<ul style="list-style-type: none">• Use to reverse transcribe RNA-containing complementary sequences only

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2) La retrotranscripción

RT Master Mix:

Step	Action																																						
1.	In a 0.2-mL microcentrifuge tube, prepare a reaction mix for all total RNA samples to be reverse transcribed. If preparing four samples, follow the recommended volumes shown below.																																						
	<table border="1"> <thead> <tr> <th rowspan="2">Component</th> <th colspan="2">Volume (μL)</th> <th rowspan="2">Final Conc.</th> </tr> <tr> <th>Per Sample</th> <th>Reaction Mix (x4)</th> </tr> </thead> <tbody> <tr> <td>RNase-free water</td> <td>Variable[‡]</td> <td>Variable[‡]</td> <td>—</td> </tr> <tr> <td>10X RT Buffer</td> <td>1.0</td> <td>4.0</td> <td>1X</td> </tr> <tr> <td>25 mM MgCl_2</td> <td>2.2</td> <td>8.8</td> <td>5.5 mM</td> </tr> <tr> <td>deoxyNTPs Mixture (2.5 mM)</td> <td>2.0</td> <td>8.0</td> <td>500 μM per dNTP</td> </tr> <tr> <td>Random Hexamers (50 μM)</td> <td>0.5</td> <td>2.0</td> <td>2.5 μM</td> </tr> <tr> <td>RNase Inhibitor (20 U/μL)</td> <td>0.2</td> <td>0.8</td> <td>0.4 U/μL</td> </tr> <tr> <td>MultiScribe Reverse Transcriptase (50 U/μL)</td> <td>0.625</td> <td>2.5</td> <td>3.125 U/μL</td> </tr> <tr> <td>Total[§]</td> <td>6.525</td> <td>26.1</td> <td>—</td> </tr> </tbody> </table> <p>‡ The volume of RNase-free water (μL) will be $3.475 - \text{RNA sample volume}$ in a 10-μL reaction.</p> <p>§ If changing the reaction volume, make sure the final proportions are consistent with the recommended values above.</p>	Component	Volume (μL)		Final Conc.	Per Sample	Reaction Mix (x4)	RNase-free water	Variable [‡]	Variable [‡]	—	10X RT Buffer	1.0	4.0	1X	25 mM MgCl_2	2.2	8.8	5.5 mM	deoxyNTPs Mixture (2.5 mM)	2.0	8.0	500 μM per dNTP	Random Hexamers (50 μM)	0.5	2.0	2.5 μM	RNase Inhibitor (20 U/ μL)	0.2	0.8	0.4 U/ μL	MultiScribe Reverse Transcriptase (50 U/ μL)	0.625	2.5	3.125 U/ μL	Total[§]	6.525	26.1	—
Component	Volume (μL)		Final Conc.																																				
	Per Sample	Reaction Mix (x4)																																					
RNase-free water	Variable [‡]	Variable [‡]	—																																				
10X RT Buffer	1.0	4.0	1X																																				
25 mM MgCl_2	2.2	8.8	5.5 mM																																				
deoxyNTPs Mixture (2.5 mM)	2.0	8.0	500 μM per dNTP																																				
Random Hexamers (50 μM)	0.5	2.0	2.5 μM																																				
RNase Inhibitor (20 U/ μL)	0.2	0.8	0.4 U/ μL																																				
MultiScribe Reverse Transcriptase (50 U/ μL)	0.625	2.5	3.125 U/ μL																																				
Total[§]	6.525	26.1	—																																				

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


Samples and endogenous controls are arranged on three plates as shown below. 50 μL of PCR master mix containing cDNA are added to each well.

The reactions are kept on ice until the plates are loaded on the 7300/7500 system.

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4) Análisis de los datos (*7500 Software v 2.0.6*).

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7500 Software v 2.0.6

QuantStudio Design & Analysis Software

En la nube de ThermoFisher:

(<https://www.thermofisher.com/es/es/home/cloud.html>)

DA Design and Analysis Application.

RQ Relative Quantification (análisis multiplaca).